

QROXB2® PERIO FotoSan Blue Agent - Instructions for use (IFU)

Product description

QROXB2® FotoSan® Blue Agent is available in two viscosities:

- QROXB2® ENDO FotoSan Blue Agent for use in root canals.
- QROXB2® PERIO FotoSan Blue Agent for use in periodontal pockets and on surfaces.

Both products contain the same active substances at the same concentration. The only difference is the viscosity. The PERIO formulation is presented as a gel to facilitate application in periodontal pockets and on surfaces.

Instructions for use (IFU)

Package contents

- 1 syringe containing Hydrogen Peroxide Gel
- 1 syringe containing QROXB2 Booster
- 1 Luer-lock connector
- 1 application tip

Preparation

1. Remove the caps from both syringes.
2. Attach the Luer-lock connector to one of the syringes.
3. Connect the two syringes securely.
4. Transfer the contents from one syringe to the other by fully depressing the plunger.
5. Mix the components by alternately transferring the material back and forth between the two syringes.
6. Repeat the mixing process **15 times**.
7. Transfer the mixed gel into the syringe intended for application.
8. Remove the Luer-lock connector.
9. Attach the application tip.

Use

The product is now ready for use according to the clinical protocol for QROXB2® PERIO FotoSan Blue Agent.

Important

Use immediately after mixing.



CMS Dental A/S, H C Ørstedes Vej 16B, 7800 Skive, Denmark

QROXB2® ENDO FotoSan Blue Agent - Instructions for use (IFU)

Product description

FotoSan® Blue Agent is available in two viscosities:

- QROXB2® ENDO FotoSan Blue Agent for use in root canals.
- QROXB2® PERIO FotoSan Blue Agent for use in periodontal pockets and on surfaces.

Both products contain the same active substances at the same concentration. The only difference is the viscosity. The PERIO formulation is presented as a gel to facilitate application in periodontal pockets and on surfaces.

Instructions for use (IFU)

Package contents

- 1 syringe containing Hydrogen Peroxide Liquid
- 1 syringe containing QROXB2 Booster
- 1 Luer-lock connector
- 1 application tip

Preparation

1. Remove the caps from both syringes.
2. Attach the Luer-lock connector to one of the syringes.
3. Connect the two syringes securely.
4. Transfer the contents from one syringe to the other by fully depressing the plunger.
5. Mix the components by alternately transferring the material back and forth between the two syringes.
6. Repeat the mixing process **10 times**.
7. Transfer the mixed solution into the syringe intended for application.
8. Remove the Luer-lock connector.
9. Attach the application tip.

Use

The product is now ready for use according to the clinical protocol for QROXB2® ENDO FotoSan Blue Agent.

Important

Use immediately after mixing.

