

User tips for mixing the ProlacSan in the syringes



1. Open the foil bage with scissors.



2. Withdraw the punger to the 1.0 ml if the plunger has not already been pulled back.



3. Flick the syringe to distribute and loosen the powder.



4. Draw up 1.5 ml water until the plunger is between 2.0 and 2.5 ml. Shake vigorously until all the powder is dissolved. This takes approx. 1 minute



5. Depress the plunger to release air. Wait for gel to form, approx. 15 minutes. If the steel ball blocks the injection of the gel, withdraw plunger, shake and continue.



6. Because of the steel ball the plunger cannot be depressed completely. To use the remaing gel suck in air and depress the plunger again.

The patient should refrain from drinking or eating for a minimum of 30 minutes. Follow up with ProlacSan tablets at home, 1 daily.